



## RAFR05204 Advanced Under Counter Freezer

| Part Number                 | RAFR05204                                   |
|-----------------------------|---|
| Capacity (Litres)           | 115   |
| Style                       | Under Counter                               |
| Exterior Dimensions (HxWxD) | 835 x 600 x 600mm                           |
| Energy Consumption          | 1.5kWh/24h                                  |
| Operating Temperature       | -20°C                                       |
| Temperature Set Point Range | -18°C to -22°C                              |
| Interior Construction       | Stainless Steel                             |
| Exterior Construction       | Stainless Steel                             |
| Shelves                     | 2   |
| Shelf dimensions (W x D)    | 460 x 4075mm<br>460 x 240mm                 |
| Door Hinging                | Right<br>(can be reversed at time of order) |
| Self Closing Door           | ✓   |
| Door Lock                   | ✓   |
| Levelling Feet              | 4   |
| Access Port                 | ✓   |
| Remote Alarm Contacts       | ✓   |
| Digital Temperature Control | ✓   |
| Digital Temperature Display | ✓   |
| High Temperature alarm      | ✓   |
| Low Temperature alarm       | ✓   |
| Door Open alarm             | ✓   |
| Manual defrost              | ✓   |
| GWP Friendly Refrigerant    | ✓   |
| Ammonia Free                | ✓   |
| 2 Year Warranty*            | ✓   |

Robustly constructed with a stainless Steel exterior this freezer from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products .



Insulated with Polyurethane Water Base Insulation the refrigerant is also ammonia, HCFC and CFC free and uses a GWP friendly refrigerant.

Fitted with a sturdy door lock, this Advanced freezer is ideal for demanding environments such as busy laboratories and features a manual defrost putting you in control.

Remote alarm contacts and an access port are also fitted as standard.

\*Parts only warranty outside of the UK

