



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Labcold Cooled Incubators

Cooled Incubator 150L Solid Door RLSD01503

LABCOLD Cooled Incubators have controlled cooling and heating allowing you to choose an operating temperature between +5°C and +40°C unlike normal incubators which can only operate from several degrees above room temperature.

Cooled incubators are usually a necessity for reliable incubation at or below 30°C in most laboratory environments.

Specifications

Product Code	RLSD01503
Gross Capacity	150 litres
Style	Undercounter
UK Warranty	1 year
External Dimensions (mm) HxWxD	835 x 600 x 600
Internal Dimensions (mm) HxWxD	690 x 480 x 410
Net Weight	43kg
Temperature Range	5°C to 40°C

Features

- Digital temperature display
- Control stability: $\pm 0.1^{\circ}\text{C}$
- Door lock
- Internal light
- Glass door model is RLCG01503
- Temperature uniformity: better than $\pm 0.5^{\circ}\text{C}$ at 4°C and $\pm 1^{\circ}\text{C}$ at 20°C
- Adjustable shelving
- Fan circulated interior
- Optional internal power socket

Models

	Model	Gross Capacity	Door Style	Lock Type
	Cooled Incubator 150L Solid Door RLSD01503	150 litres	Solid	Key lock
	Cooled Incubator Glass Door 150L RLCG01503	150 litres	Glass	Key lock